

DOWNLOAD

Determination of the (15n14n) of Nitrate in Solids: Rsil Lab Code 2894

By -

BiblioGov. Paperback. Book Condition: New. This item is printed on demand. Paperback. 38 pages. Dimensions: 9.7in. x 7.4in. x 0.1in.The purpose of Reston Stable Isotope Laboratory (RSIL) lab code 2894 is to determine the (15N14N), abbreviated as 15N, of nitrate (NO3-)in solids. The nitrate fraction of the nitrogen species is dissolved by water (called leaching) and can be analyzed by the bacterial method covered in RSIL lab code 2899. After leaching, the 15N of dissolved NO3- is analyzed by conversion of NO3- to nitrous oxide (N2O), which serves as the analyte for mass spectrometry. A culture of denitrifying bacteria is used in the enzymatic conversion of NO3- to N2O, which follows the pathway shown in equation 1. Because the bacteria Pseudomonas aureofaciens lacks N2O reductive activity, the reaction stops at N2O, unlike the typical denitrification reaction, which goes to N2. After several hours, the conversion is complete, and the N2O is extracted from the vial, separated from water vapor by Nafion drier and from CO2 with a layered Mg(ClO4)2Ascarite trap, and trapped in a small-volume trap immersed in liquid nitrogen. After the N2O is released, it is further purified by gas chromatography (GC) before introduction to the isotope-ratio mass spectrometer...



Reviews

The ideal ebook i possibly study. Better then never, though i am quite late in start reading this one. It is extremely difficult to leave it before concluding, once you begin to read the book.

-- Ava Witting

The ideal ebook i possibly study. Better then never, though i am quite late in start reading this one. It is extremely difficult to leave it before concluding, once you begin to read the book.

-- Ava Witting